

Battery powered

Portable Generator



1040BT



90BT Plus



90BT



9020BT



1025BT

Flat Panel Detector with acquisition software



mini90



Scintillator: CsI
 Active array: 2560 * 3072 pixels
 Active area (H * V): 358.4mm
 Pixel pitch: 0.14 (140 μm)
 Weight: 2.9kg (including battery pack)

*DR acquisition SW with DICOM compatibility

Technical Specifications

Model		mini90	90BT	9020BT	90BT Plus	1025BT	1040BT
Output	kV Range	40-90kV				40-100kV	
	mA Range	8-10mA	12-20mA		12-30mA	15-25mA	10-40mA
Tube	Focal Spot	0.8 x 0.8mm	1.2 x 1.2mm	0.8 x 0.8mm		1.2 x 1.2mm	
Collimator	Lamp	LED					
Power	Input Power	110V-240VAC					
APR		16	12 or 7 mode				
Net Weight		4kg	8.5kg	9.8kg		11.4kg	13.5kg

AC powered

Portable Generator



1060HF

1040HF

9020HF



Medical



VET



ULTRA 100 (5kW)

Model		9020HF	1040HF	1060HF	ULTRA100
Output	kV Range	40-90kV	40-100kV		40-110kV
	mA Range	12-20mA	10-40mA	25-60mA	20-100mA
Tube	Focal Spot	1.2 x 1.2mm		2.0 x 2.0mm	1.4 x 1.4mm
Collimator	Lamp	Halogen			LED
	Type	Slide	Knob		
Power	Input Power	110V-240VAC		220-240VAC	110-240VAC
Features	APR	24			96
	Dimensions (mm)	280*140*170	355*191*180	438*202*191	438*205*208
	Net Weight	7.5kg	11.5kg	12.9kg	14kg