

# Single Panel DR

A straight arm DR equipped with one fixed type 17 X 17 FPD provides high quality images in order to meet general and specialized diagnostic needs. Unique designed straight arm enables all types of diagnostic procedures from chest to table examination in small practices or primary / urgent care institutions .

- Motorized movement for SID control and Arm elevating (up/down).
- Easy to maintain Chest examinations – covering entire wall stand applications
- Provide finest image resolution through Flat Panel TFT detector
- The collapsible column can shrink upto 1.4m according to the space condition.
- Easy installation and operation in a small clinic or a vehicle



Detector angle adjustable




Easy Lock-free mechanics



[ “ECOVIEW 9 / Straight Arm” is a complete solution for Chiropractic / Urgent & Primary Care. ]



Table position



Axial projection



Quick positioning to Chest & Table